


Universal Base

Workflow chart



Abutments, made your way



Quickly create quality abutments for Nobel Biocare implants using whichever workflow you prefer. The whole system – from the implant and the clinical screw to the universal base – are tested together. This ensures all the components you receive, work together for optimal long-term performance.

Biomechanical investigations and micro gap measurements with cross-sectional SEM images highlight the precise fit that's only guaranteed by the original abutment-implant interface.

Selecting an abutment with a precise fit is decisive for system performance, as this ensures that occlusal forces are distributed evenly and that uncontrolled peak stresses are avoided.

Conventional press workflow

What do I need?




Clinical screw and burn-out coping included

38213	Universal Base CC NP 1.5mm
38214	Universal Base CC RP 1.5mm
38215	Universal Base CC WP 1.5mm
38216	Universal Base CC NP 3mm
38217	Universal Base CC RP 3mm
38218	Universal Base CC WP 3mm

Order online:
store.nobelbiocare.com


Nobel Biocare Universal Base and burn-out coping

How does it work?



1 Universal Base preparation

- Hand-tighten the Universal Base onto the model. Make sure to use the laboratory screw.




2 Burn-out coping preparation

- Seat the burn-out coping onto the Universal Base.
- Adjust the height of the burn-out coping according to the required occlusal plane. Make sure the Universal Base remains fully covered.



3 Production

- Create a wax-up restoration and use the standard procedure to either press or cast the coping or full-contour crown. Finalize the restoration.




4 Finalization

- Seal the screw channel with wax.
 - Sandblast the contact surface of the Universal Base with aluminium oxide 50µm at a maximum of 2 bar.
- Caution:** Do not sandblast the seating area.
- Tip:** during the blasting procedure, use a protection analog to prevent any modification of the abutment / implant interface.
- Clean the surface with alcohol or steam.
 - Bond the restoration to the Universal Base according to the cement manufacturer's instructions. Use only self-adhesive dental cement/bonding material suitable for zirconium dioxide ceramics or PMMA.

CAD/CAM workflow with exocad® software

What do I need?



Any CAD/CAM system that integrates with exocad® software



Order online:
store.nobelbiocare.com

NobelProcera Abutment Position Locator
Conical Connection NP 36736
Conical Connection RP 36737
Conical Connection WP 37569



Download here:
nobelbiocare.com/stl


Nobel Biocare stl file exocad



Order online:
store.nobelbiocare.com

Nobel Biocare Universal Base


How does it work?



1

Scanning

- Insert the NobelProcera Abutment Position Locator onto the model.
- Scan by following the usual scan routines.



2

Designing


- In the implant library select the height and implant diameter for the Universal Base.
- Design the implant crown using standard CAD tools.



3

Milling

- Send the design file to a milling unit or production facility that accepts exocad design files.



4

Finalization

- Once the restoration is produced, finalize it following the material manufacturer's instructions.
- Seal the screw channel with wax.
- Sandblast the contact surface of the Universal Base with aluminium oxide 50 µm at a maximum of 2 bar.
Caution: Do not sandblast the seating area.
- Clean the surface with alcohol or steam.
- Bond the restoration to the Universal Base according to the cement manufacturer's instructions. Use only self-adhesive dental cement/bonding material suitable for zirconium dioxide ceramics or PMMA.

CAD/CAM workflow with 3Shape® Dental System™

What do I need?



3Shape Dental System™*



Order online:
store.nobelbiocare.com

ELOS Desktop Scan Bodies

Conical connection NP	2B-A
Conical connection RP	2B-B
Conical connection WP	2B-C



Download here:
nobelbiocare.com/stl

Nobel Biocare DME file for 3Shape®



Order online:
store.nobelbiocare.com


Nobel Biocare Universal Base

How does it work?



1 Scanning

- Use the ELOS Pinol Screwdriver to screw the ELOS Desktop Scan Body on the the model.
- Scan following the usual scan routines.



2 Designing


- In the implant library, select the height and implant diameter for the Universal Base.
- Design the implant crown using standard CAD tools.



3

Milling

- Send the design file to a milling unit or production facility accepting 3Shape® design files.



4

Finalization

- Once the restoration is produced, finalize it following the material manufacturer's instructions.
- Seal the screw channel with wax.
- Sandblast the contact surface of the Universal Base with aluminium oxide 50 µm at a maximum of 2 bar.
Caution: Do not sandblast the seating area.
- Clean the surface with alcohol or steam.
- Bond the restoration to the Universal Base according to the cement manufacturer's instructions. Use only self-adhesive dental cement/bonding material suitable for zirconium dioxide ceramics or PMMA.

* To the date of June 2015. Supported 3Shape® Scanners*: D640, D700, D710, D800, D810, D900 and D900L. Supported 3Shape® Software*: 2012 (2.7.8.13), 2013 (2.8.8.4-2.8.8.9), 2014 (2.9.9.0-2.9.9.5) and 2015 (2.15.2.0-2.15.3.1).


Scanner types and software versions are continuously tested and validated. Refer to your Nobel Biocare sales office for currently supported scanner and software.

CAD/CAM workflow without component integration

What do I need?




Any CAD/CAM System, including intraoral scanning systems



Nobel Biocare Universal Base

How does it work?



1

Scanning


- Block out screw access hole and undercuts of the mounted Universal Base (e.g. with putty), in case the software does not support virtual blocking out.
- Scan spray may be applied.



2

Designing


- Design the implant crown using standard CAD tools.



3

Milling

- Send the design file to a milling unit or production facility that accepts the design files.




4

Finalization

- Once the restoration is produced, finalize it following the material manufacturer's instructions.
- Seal the screw channel with wax.
- Sandblast the contact surface of the Universal Base with aluminium oxide 50 µm at a maximum of 2 bar.
Caution: Do not sandblast the seating area.
- Clean the surface with alcohol or steam.
- Bond the restoration to the Universal Base according to the cement manufacturer's instructions. Use only self-adhesive dental cement/bonding material suitable for zirconium dioxide ceramics or PMMA.

Clinical procedure

How does it work?



1

Preparation

- Remove the restoration from the implant analog.
- Clean and disinfect according to the cleaning and sterilization guidelines.*

* refer to nobelbiocare.com/sterilization

2

Insertion

- After removing the temporary abutment or healing cap, connect and tighten the final restoration to the implant.
- It is recommended to verify the seating of the restoration with a radiographic scan.
- Tighten the Universal Base restoration to 35Ncm using the Unigrip Screwdriver and Manual Torque Wrench Prosthetic.

Important: Never exceed 35Ncm prosthetic tightening torque for the abutment screw. Overtightening of the abutment may lead to screw fracture.


3

Finalization

- Block-out the screw head (e.g. with teflon, gutta-percha, cotton, any permanent soft resin etc.) before closing the screw access hole with composite. This ensures that no composite blocks access to the screw head, allowing easy removal of the restoration at any time.
- Close the screw access hole using dental composite, adhering to the manufacturer's bonding and curing guidelines.



Be optimally equipped with the prosthetic kit



Prosthetic kit 37448