

# On1™ concept

## Product overview

**On1 is a new restorative concept. It preserves the soft tissue attachment by acting as a base for abutments from the time of implant placement to finalization. This leaves the soft tissue undisturbed during impression taking and placement of restorative components.**

### Restore implants with ease

Connecting the impression coping or your preferred abutment is simple, as the base places the connection at tissue level.

### Choose between cement or screw-retention


Depending on indication or clinical preference, you can choose a cement- or screw-retained final restoration.

### Preserve the soft tissue attachment


Leave the soft tissue undisturbed for optimized healing. The base remains in position while the healing, temporary and final abutments are placed.

### Retain restorative flexibility


The two base-height options make it possible to optimize short- and long-term esthetic outcomes.



On1 Base




Healing



Temp/Final


# One base, many restorative options

## Surgical components


On1 Base incl. handle				
		NP	RP	WP
	1.75 mm	38689	38691	38693
	2.5 mm	38690	38692	38694

Delivered with the corresponding On1 Clinical Screw


STERILE

On1 Healing Cap				
		NP	RP	WP
	1.5 mm	38695	38697	38699
	2.5 mm	38696	38698	38700


STERILE

On1 IOS Healing Cap incl. handle				
		NP	RP	WP
	4 mm			38739
	4.5 mm	38735	38737	
	5 mm			38738
	6 mm	38734	38736	

STERILE


On1 Clinical Screw				
		NP	RP	WP
	1.75 mm	38744	38746	38748
	2.5 mm	38745	38747	38749

STERILE


On1 Base Screwdriver		
	Machine 21 mm	38752
	Manual 25 mm	38753




## Lab components

On1 Base Replica			
	NP	RP	WP
	38731	38732	38733




On1 IOS Replica			
	NP	RP	WP
	38740	38741	38742




On1 Prosthetic Lab Screw			
	NP	RP	WP
	38751	38751	38751



On1 IOS Insertion Tool			
	NP	RP	WP
	38201	38201	38201



## Easy handling



The On1 Base, IOS Healing Cap and Temporary Abutment are all delivered with a pre-mounted handle for easier placement.

## Restorative components

On1 Temporary Abutment incl. handle			
	NP	RP	WP
<b>Engaging</b>	38701	38703	38705
<b>Non-engaging</b>	38702	38704	38706

Delivered with the corresponding On1 Prosthetic Screw

STERILE R

On1 Impression Coping			
	NP	RP	WP
<b>Closed Tray</b>	38725	38726	38727
<b>Open Tray</b>	38728	38729	38730



On1 Final Abutment				NP	RP	WP
<b>Universal Abutment (screw-retained solution)</b>	0.3mm	38707	38709	38711		
	1.25mm	38708	38710	38712		
<b>Esthetic Abutment Titanium</b>	0.3mm	38713	38714	38715		
<b>Esthetic Abutment Zirconia</b>	0.3mm	38716	38717	38718		

Delivered with the corresponding On1 Prosthetic Screw



On1 Prosthetic Screw			
	NP	RP	WP
	38750	38750	38750

STERILE R

Unigrip™ Screwdriver				
<b>Manual</b>	<b>20 mm</b>	<b>28 mm</b>	<b>36 mm</b>	
	29148	29149	29150	
<b>Machine</b>	<b>20 mm</b>	<b>25 mm</b>	<b>30 mm</b>	<b>35 mm</b>
	29151	29152	29153	29154



Manual Torque Wrench Prosthetic	
29165	




Manual Torque Wrench Adapter	
29167	





For the full assortment of Nobel Biocare products, see our current product catalog.

[nobelbiocare.com/catalog](http://nobelbiocare.com/catalog)



### Order online

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

[nobelbiocare.com/store](http://nobelbiocare.com/store)



### Order by phone

Call our customer service team, or contact your sales representative.



On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit [nobelbiocare.com/warranty](http://nobelbiocare.com/warranty)



[nobelbiocare.com/on1](http://nobelbiocare.com/on1)

80140 GB 1604 Printed in the EU © Nobel Biocare Services AG, 2016. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](http://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. Disclaimer: Some products may not be regulatory cleared/ released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed dentist. See Instructions for Use for full prescribing information, including indications, contraindications, warnings and precautions.